

Oxygen, anytime, anywhere.

INOGEN ONE^{G2}



Oxygen therapy should not limit you in any way.

With the INOGEN ONE G2 there are no worries about running out of oxygen, no need to keep an eye on the clock and no need to carry around heavy equipment. Simply use the INOGEN ONE G2 wherever you are, wherever you go and whenever you want.

Restore your freedom with the INOGEN ONE^{G2}

INOGEN ONE G2 Whenever you want
INOGEN ONE G2 is built to withstand long term and heavy duty use. You can expect it to run for years and it will only need full maintenance after it has been running for 30,000 hrs.

INOGEN ONE G2 Wherever you go
Restoring your freedom is definitely a large part of what makes the INOGEN ONE G2 so special, so go ahead, jump in the car to run errands or take a weekend trip to see your family, or take that long trip you never dreamed possible. Head home when you want, not when lack of oxygen forces you to.

INOGEN ONE G2 Wherever you are
It sits discreetly at your feet, behind a chair or next to the bed and it's light enough to carry over your shoulder. All of your oxygen needs can be met with this one small device, day or night, home or away. When using the INOGEN ONE G2 you will embrace the benefits of this small and lightweight device, producing just a very light sound, hardly noticeable to you or those around you. And, it's easy to operate too, just a couple of buttons.

Most airlines allow you to use the INOGEN ONE G2 when you fly, making it easier than ever to reclaim your independence.

Technical Specifications:

Weight :	3.2 kg with battery.
Size:	Height: 24cm; Width: 10cm; Length: 27cm.
Oxygen Flow:	Pulse dose delivery system. Five flow settings from 1-5.
Power Sensing:	Universal Power Supply with auto-sensing for worldwide use. In AC mode for 100- 240V; in DC mode for 12V for use in the car/airplane/ boat/camper.
Battery:	Duration 3.5-4 hours*, Optional 7-8 hour battery.
Noise Level:	37 Decibels.*
Warranty:	3 year limited warranty.**
Operation:	Simple control functions; clear LCD display.
Use:	Stationary, portable and travel. Clinically effective for day and night use.

*At setting 2, in typical home use.

Accessories Included



Lithium Ion Battery

The battery lasts for up to four hours. When it's plugged in (AC or DC), the battery automatically recharges. Eight hour batteries are also available.



Universal Power Supply

When plugged into an AC or DC power outlet, the INOGEN ONE G2 delivers an unlimited supply of oxygen.



Cart & Carry Bag

Designed specifically for the INOGEN ONE G2. Take advantage of your new independence and mobility.



sales@smilerehab.com; www.oxygenconcentrators.co.uk

Ability House, 109 New Greenham Park, Newbury, Berkshire, RG19 6HN

Call us on
01635 37550